

**THERAPEUTIC CAP-
BONS** must be used
under a physician's
direction.

**CAUTION: PROTECT
EYES** with dark glass,
light-tight goggles.

EXCESS EXPOSURE
may cause severe sun-
burn.

EXAMINE package and stamping on each carbon to be sure that prescription is properly filled. Types are identified by letters, sizes by number.

"B", "C" and "K" carbons radiate intense ultra-violet light.

"E" carbons radiate infra-red light, and less ultra-violet light.

Indicated initial exposure for low amperage arcs is not more than three (3) minutes at not less than three (3) feet from the lamp. For arcs of higher amperage, or in case of unusual susceptibility to sunburn, duration should be proportionately less. Increase subsequent exposures gradually.

PRINTED IN U. S. A.

"B," "C" and "K" CARBONS radiate intense ultra-violet light.

"E" CARBONS radiate infra-red and less ultra-violet light.

READ INSTRUCTIONS on INSERT in package.

THERAPEUTIC CAR-
BONS must be used
under a physician's
direction.

CAUTION: PROTECT
EYES with dark glass,
light-tight goggles.

EXCESS EXPOSURE
may cause severe sun-
burn.

EXAMINE package and stamping on each carbon to be sure that prescription is properly filled. Types are identified by letters, sizes by number.

"B", "C" and "K" carbons radiate intense ultra-violet light.

"W" carbons radiate infra-red light, and less ultra-violet light.

Indicated initial exposure for low amperage arcs is not more than three (3) minutes at not less than three (2) feet from the lamp. For arcs of higher amperage, or in case of unusual susceptibility to sunburn, duration should be proportionately less. Increase subsequent exposures gradually.

PRINTED IN U. S. A.

Use SUPER-TAN Carbons (marked with a RED line on the butt end) for an initial exposure of not more than three (3) minutes at not less than three (3) feet from the lamp.

Use SUNSHINE Carbons by starting at five minutes at three feet.

Increase subsequent exposures gradually. Do not decrease distance from lamp.

These carbons are not intended for therapeutic use. These suggestions are for persons of average susceptibility to sunlight and under normal conditions. Under other conditions CONSULT YOUR PHYSICIAN.

PRINTED IN U. S. A.

CAUTION: PROTECT EYES with dark glass, light-tight goggles.

QUANTITY	TYPE	SIZE
50	C	3 mm. x 6 in.

H. G. FISCHER & CO., Inc.

2323-2337 WABANSA AVENUE
CHICAGO, ILLINOIS

1774
SUNSHINE CARBONS

SUNSHINE CARBONS

THERAPEUTIC CARBONS

THERAPEUTIC CARBONS



TRADE MARK

© 1928 N. C. CO. INC.

SUPER-TAN CARBONS

SUPER-TAN CARBONS

DIRECTIONS

FOR USERS of SUNSHINE and SUPER-TAN CARBONS
PROTECT YOUR EYES; WEAR DARK GLASS, LIGHT-TIGHT GOGGLES
TIME YOUR SUN-BATH TO AVOID SEVERE SUNBURN.

IBI Sanamore 医療用カーボン
THERAPEUTIC CARBONS D

(株) 東京光線療法研究所 8×150mm
10本入

IBI Sanamore 医療用カーボン
THERAPEUTIC CARBONS B

(株)東京光線療法研究所 8×150mm
10本入

DIRECTIONS

CAUTION:

THERAPEUTIC CARBONS MUST *be* USED *as* DIRECTED *by a* PHYSICIAN.

PERSONS EXPOSED SHOULD WEAR DARK GLASS, LIGHT-TIGHT GOGGLES.

EXCESS EXPOSURE MAY CAUSE SEVERE SUNBURN.

DIRECTIONS

CAUTION:

THERAPEUTIC CARBONS MUST *be* USED *as* DIRECTED *by a* PHYSICIAN.

PERSONS EXPOSED SHOULD WEAR DARK GLASS, LIGHT-TIGHT GOGGLES.

EXCESS EXPOSURE MAY CAUSE SEVERE SUNBURN.

DIRECTIONS

CAUTION:

FOR USERS *of* SUNSHINE *and* SUPER-TAN CARBONS

PROTECT YOUR EYES, WEAR DARK GLASS, LIGHT-TIGHT GOGGLES.

TIME YOUR SUN-BATH TO AVOID SEVERE SUNBURN.

DIRECTIONS

Use SUPER-TAN Carbons (marked with a RED line on the butt end) for an initial exposure of not more than three (3) minutes at not less than three (3) feet from the lamp.

Use SUNSHINE Carbons by starting at five minutes at three feet.

Increase subsequent exposures gradually. Do not decrease distance from lamp.

These carbons are not intended for therapeutic use. These suggestions are for persons of average susceptibility to sunlight and under normal conditions. Under other conditions CONSULT YOUR PHYSICIAN.

PRINTED IN U. S. A.

A Daily Sun-Bath Promotes Health

MADE IN U.S.A. NATIONAL CARBON CO., INC. CLEVELAND, OHIO
Sole U.S. Agent
Cleveland, Ohio and Carbon Corporation, San Francisco, California

For a Quick Healthy Tan

Designed for Use in Light Therapy

A Daily Sun-Bath is Good for You

Intended for Physicians' Use

For a Quick Tan

THERAPEUTIC CARBONS

QUANTITY

50

TYPE

C

SIZE

8 mm. x 6 in.

H. G. FISCHER & CO., Inc.

2323-2327 WABANSIA AVENUE

CHICAGO, ILLINOIS

"BLUE SUN (LOEB) DUO CORE"

SUPER-TAN

THERAPEUTIC

"YELLOW DEIOXICANI (LOEB) DUO CORE"

ACTINO LABORATORIES, INC.

"YELLOW SUN (LOEB) DUO CORE"

ACTINO LABORATORIES, INC.

"YUNGSTEN (LOEB) DUO CORE"

ACTINO LABORATORIES, INC.

PRISMATIC

"IN (LOEB) DUO CORE"

MICRO PROJECTOR

PRISMATIC

SPECTRUM (LOEB)

ACTINO LABORATORIES, INC.

ACTINO LABORATORIES, INC.

THERAPEUTIC

THERAPEUTIC

THERAPEUTIC

THERAPEUTIC

SPECTRUM (LOEB)

ACTINO LABORATORIES, INC.

